

Fluorescence Fingerprinting Protocol (Tara)

Pw: eclipse

➔ Scan

➔ Setup

- Cary – Em. Excitation 200 3D mode: Check (450 & 10)
Start: 250 Medium
Stop: 600
- Options – PMT = high
- Accessories – Multicell #2
- Autostore – On: Prompt on Start (ASCII csv)
- OK
- Start with blank: MilliQ and press Start
 - Only want to see Raleigh Peaks
- When scan is complete click on graph then box at top to see 3D plot

UV-Vis

Pw: cary50

* corrects for inner filtering effects for PARAFAC (inner filtering corrections required if $A > 0.3$ at 254 nm)

➔ Scan ➔ Connect

➔ Setup

- Start: 600 Stop: 200 Scan Control: Medium
 - Baseline: zero/baseline correction
 - Reports: Select for ASCII (csv)
 - Autostore: On: prompt on start
 - Baseline (to start scan)
- ➔ On prompt add onecard to block beam, once done remove card and press START

Cleaning Cuvette:

1. Rinse with MilliQ and pat dry
2. Rinse with Standard/Sample a couple times then fill

Standard:

10 mg C/L Southampton DOC

2.4 μ M Tyrosine

1 μ M Tryptophan

Cuvette Preparation and Storage:

- Store in MilliQ plus cuvette soap
- Before use put in beaker with MilliQ then add soap
- Heat for at least 10 min in oven @ 40°C

Data Analysis:

- Save data from instrument computers to USB
- Remove text portion (rows) and wavelength (columns)
- 1st box arbitrary value of -99; 2nd box start @ 200 and increase by 10 across row (should end at 450)
- Resave (as short a name as possible)
- Move to same folder as matlab scripts to run analysis
- 10 points to remove Raleigh scattering